pegAssist™
Insole

USA/GB Indication
<ul style="list-style-type: none"> > Selective off-loading > Acute treatment > Pain therapy > Trauma care > Postoperative care
Contraindication
<ul style="list-style-type: none"> > Open fractures
Technical data
<p>Size:</p> <ul style="list-style-type: none"> > Seven sizes available > Right and left fit <p>All required sizes available for the following foot and ankle products:</p> <ul style="list-style-type: none"> > OrthoWedge – Off-loading shoe > OrthoWedge Light – Off-loading shoe > HeelWedge® – Off-loading shoe > MedSurg™ – Healing shoe > APB® – All purpose boot > SlimLine® – Cast boot > FXPro® Walker – Stabilising orthosis <p>Material: The 18 mm thick insole consists of 3 layers: a combination of EVA, Plastizote® and a top layer made of Poron®.</p> <p>Maintenance: The product requires no maintenance.</p>
Cleaning instructions
<p>The soles can be cleaned with warm water and soap or regular disinfectants.</p>

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PegAssist™ Insole

Off-loading Insoles

The PegAssist™ Insole is a system effectively off-loading the plantar aspect of the foot. The easy removing of the single pegs allows an individual customization in the phase of postoperative care or acute treatment. The thin cardboard cover stabilizes the off-loaded area and prevents the surrounding pegs from shifting. The pressure is reduced up to 60%.

Recommendations for use:

1. Remove the black original insole from the shoe (fig. 1).
2. Insert in return the PegAssist™ Insole into the shoe with the blue side facing upwards (fig. 2).
3. Encircle the wound or pressure ulcer with the enclosed marker pen on the sole of the patient (fig. 3).
4. The patient should now stand up in the shoe onto the insole which then pictures the marked area onto the blue upper layer of the insole (fig. 4).
5. Carefully remove the black pegs beneath the marked area where the off-loading is taking place (fig. 5).
6. Apply the complete surface of the cardboard cover (included in delivery) onto the black insole. The cover provides an adhesive backing which affixes the pegs surrounding the off-loaded area and ensures non shifting (fig. 6).
7. Reinsert the adjusted PegAssist™ Insole into the shoe (fig. 7).
8. Ready (fig. 8).

Caution

- > It is recommended to wear bandages or socks.

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