# Body Armor® Ankle Flex

- EN Ankle Bandage
- DE Sprunggelenkbandage
- ES Tobillera
- FR Cheville
- IT Cavigliera
- NL Teenspalk
- DA Hallux valgus-tâskinne
- SV Fotledsförband
- PL Bandaż stawu skokowego
- SL Povoij za gleženj
- HR Zavoj za skočni zglob
- PT Bandagem para tornozelo

# DARCO

## Body Armor® Ankle Flex

1.

2.

## EN Instructions for use

### Body Armor® Ankle Flex | Ankle Bandage

#### Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This aid assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product on injured, irritated, or damaged skin. This aid is intended for single-patient use only.

## Purpose

The Body Armor® Ankle Flex stabilises and protects the ankle according to the indication.

## Indications / Contraindications

### Indications

- > Osteoarthritis of the upper and lower ankle joint
- > Rheumatic diseases
- > Prophylactic use during daily life for ligament weakness
- > Postoperative and post-traumatic irritation
- > Articular effusion and swelling caused by osteoarthritis or arthritis
- > Sense of instability

### Contraindications

- > Skin changes due to allergic reactions, inflammation or injury (e.g. swelling, redness)
- > Impaired circulation or lymphatic soft tissue swelling.
- > Neurogenic impairments of the sensory function and trophic condition of the skin in the area of the body to be treated (sensory impairment with and without skin damage)
- > Impaired perfusion and sensation in the foot, e.g. diabetes mellitus

## Design features

Areas of different strength, silicone pads incorporated in the medial and lateral regions and plantar reinforcement support the foot's physiological structures.

## EN Use

### Putting on the bandage

- > The opening to which the DARCO flag is attached is at the top.
- > Roll the upper part of the bandage down to the upper edge of the silicone inserts.
- > Pull the bandage over the foot until the heel is properly positioned in the anatomically knitted heel part. (1)
- > The silicone inserts should then be positioned behind and below the medial and lateral malleolus.
- > Roll the upper part of the bandage back up and pull the bandage into the correct position. (2)
- > The Body Armor® Ankle Flex is properly positioned when the seam is centred over the sole of the foot and the heel.

## Combination with other products

No specifications


## Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling/cleaning.

## Technical data

### Sizes:

	Ankle measurement in cm	inches
XS	< 19	< 7 <sup>1</sup> / <sub>2</sub>
S	19 – 21	7 <sup>1</sup> / <sub>2</sub> – 8 <sup>1</sup> / <sub>4</sub>
M	21 – 23	8 <sup>1</sup> / <sub>4</sub> – 9
L	23 – 25	9 – 9 <sup>7</sup> / <sub>8</sub>
XL	25 – 27	9 <sup>7</sup> / <sub>8</sub> – 10 <sup>7</sup> / <sub>8</sub>

Can be worn on the right or left foot

**Colour:** Grey-Blue

**Material:** 65% polyamide, 24% elastodiene/latex, 8% viscose, 3% elastane

## Maintenance

The product does not require maintenance.

## Cleaning instructions

- > Hand wash (cold or lukewarm, max. 40°C)
- > Do not bleach
- > Do not dry in dryer
- > Do not iron
- > Do not dry-clean

## Disposal

Dispose of after use Please dispose of the packaging in an environmentally manner.  
The product can be disposed of with household waste.

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