# MECRON Hip-Spica

- EN Hip Spica Bandage
- DE Hüftspicaverband
- ES Vendaje en espica para cadera
- FR Bandage de compression de la hanche
- IT Bendaggio a spiga per anca
- NL Spicaverband voor de heup
- DA Hoftespica
- SV Höftförband med 8-turer
- PL Opatrunek klasowy biodra
- SL Kolčna klasasta obveza
- HR Spica zavoj za kuk
- PT Ligadura de imobilização da anca

**DARCO**

**MECRON**  
*by* **DARCO**

## EN Instructions for use

### MECRON Hip-Spica | Hip Spica Bandage

#### Important notes

If there is a serious deterioration in your state of health as a result of using the medical device, please report this to your specialist dealer or the manufacturer as soon as possible. If you report this to the specialist dealer, the dealer will forward this immediately. The reports are agreed with and processed by the competent authorities.

Please always follow these instructions for use. In case of adverse effects (e.g. skin swelling, pressure points, disorders of blood circulation etc.) please tell your doctor. Function and optimum protection are no longer guaranteed if the product is handled incorrectly. This aid assists recovery and should be used exclusively under the instruction of and in consultation with your doctor or a healthcare professional. DARCO (Europe) GmbH does not assume any responsibility. Do not wear the product with injured, irritated, damaged or intact skin. You should use bandages, for example.

We recommend that you do not drive a vehicle while you are wearing this aid. This aid is intended for single-patient use only.

## Purpose

The hip spica bandage is used for compressing and stabilising the hip joint.

## Indications / Contraindications

#### Indications

- > Stabilization of the hip joint following total hip replacement
- > Repositioning of a dislocated hip joint
- > Hip injury
- > Chronic hip instabilities
- > Post-traumatic post-operative immobilization

#### Contraindications

No contraindications are known.

## Design features

The bandage parts are made of elastic webbing and a variably adjustable hook and loop fastener to adjust the width. The pelvic and thigh bandages are connected elastically and compress the hip area.

## Working life / Lifespan of the product

The lifespan of the medical device is determined by natural wear and tear due to proper handling / cleaning. The medical device is designed to be used for the usual duration of use prescribed by a doctor without showing a reduction in safety or performance. After the period of use prescribed by your doctor, check the function of the product together with your doctor or specialist dealer.

1. Open the hook and loop straps of the hip spica bandage.
2. Wrap the bandage around the patient's pelvis and thigh so that it fits snugly.
3. First close the hook and loop fastener on the pelvis.
4. Now you can close the hook and loop fastener on the thigh.

**Note**

Attach the bandage in such a way to ensure a snug fit. Make sure that the skin does not become trapped in the region where the hook-and-loop fasteners are. The latter should not be fastened too tightly to prevent stasis.

## Technical data

### Sizes

The hip spica bandage comes in 2 sizes depending upon circumference. Measure the hip circumference. The bandage comes in 2 versions for the left and right side and also double-sided.

Size	Circumference
M	95 – 115 cm
L	105 – 125 cm

### Material

<b>Elastic bandage:</b>	Cotton, elastane
<b>Abdominal portion:</b>	Polyurethane, polyester
<b>Yarn:</b>	Polyester
<b>Edging strip:</b>	Polyester
<b>Bungee closure system:</b>	Polyamide

**Colour:** White

## Combination with other products

None available.

## Maintenance

The product does not require maintenance.

## Cleaning instructions

The MECRON Hip Spica Bandage is washable at 30°C. The hook and loop fasteners must be closed before washing.

## Disposal

Dispose of after use. Please dispose of the packaging in an environmentally friendly manner. The product can be disposed of with household waste.

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