OPTIMA / SBi Insole 3 x 3 Heel
Modulare Innensohle

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

USA/GB Instructions for use

OPTIMA / SBi Insole 3 x 3 Heel
Modular Insole


The modular inner sole was developed to redistribute plantar pressure and is used specifically with pressure ulcers on the fore-, midfoot or heel.

- The OPTIMA / SBi Insole 3 x 3 sole system comprises nine different parts, which are arranged like a puzzle to form a three-layer inlay. For the OPTIMA Heel it is important to insert the enclosed white insole into the shoe placing it underneath the tricolor OPTIMA / SBi Insole 3 x 3 Heel.
- The soles consist of three different shore hardnesses (**red = 25 soft**, **beige = medium** and **blue = hard**) and can be combined vertically and horizontally. Simply exchanging the modules enables patients to distribute pressure optimally subject to indication. In urgent cases this can be done on an outpatient basis in a few simple steps.
- To achieve relief in an ulcerated area of the sole, you can adapt the sole elements individually.

Ulcer marking and cutting the soles:
Patient lighter than 90 kg (Mark the pressure ulcers on the beige sole!)

- Cover the pressure ulcer with a transparent plaster and use the supplied marker pen to mark the pressure ulcer on the transparent plaster on the patient's foot (Fig. 1).
- Insert the whole OPTIMA / SBi Insole 3 x 3 sole system into the shoe and ensure that the **beige** sole lies on top. (**beige – red – blue**)
- The patient should now stand on the inner sole in the shoe and walk a few paces with the shoe closed to mark the ulcer on the uppermost (**beige**) layer of the inner sole (Fig. 2). Note: Assist the patient if required.
- Now unfasten the shoe and remove the entire sole system.
- Proceed to mark the cut-out to be removed on the top side of the **beige** sole and allow about 6/8 mm of clearance to the colour marking (Fig. 3).
- Perforate the **beige** sole along the marked area with a punch tool (Fig. 4).
- Turn the **beige** sole over and use a sharp knife or the supplied scalpel to cut out the relief zone at a 45° slant along the perforation line on the unmarked side (Figs. 5 and 6).
- Now reassemble the individual sole parts according to the pressure redistribution (see diagram) and insert the individually adapted sole system into the shoe.
- Ensure that the sole you have cut forms the middle layer of the sole!
- After putting on a bandage or a sock, the shoe can now be fitted to the patient (Fig. 7).
- The hook and loop fasteners can be shortened as required.


Typical application: patient lighter than 90 kg

Patient heavier than 90 kg (Mark the pressure ulcers on the blue sole!)

- Cover the pressure ulcer with a transparent plaster and use the supplied marker pen to mark the pressure ulcer on the transparent plaster on the patient's foot (Fig. 1).
- Insert the whole OPTIMA / SBi Insole 3 x 3 sole system into the shoe and ensure that the **blue** sole lies on top. (**blue – beige – red**)
- The patient should now stand in the shoe on the **blue** inner sole and walk a few paces with the shoe closed to mark the ulcer on the uppermost (**blue**) layer of the inner sole (like in Fig. 2, but with **blue** sole). Note: Assist the patient if required.
- Now unfasten the shoe and remove the entire sole system.
- Proceed to mark the cut-out to be removed on the top side of the **blue** sole and allow about 6/8 mm of clearance to the colour marking (like in Fig. 3, but with **blue** sole).
- Place the **beige** sole beneath the **blue** sole and perforate both soles along the marked area with a punch tool (like in Fig. 4).
- Turn the **blue** sole over and use a sharp knife or the supplied scalpel to cut out Relief Zone 1 at a 45° slant along the perforation line on the unmarked side (like in Figs. 5 and 6).
- Cut out Relief Zone 2 on the **beige** sole in the same way.
- Now reassemble the individual sole parts according to the pressure redistribution (see diagram) and insert the individually adapted sole system into the shoe.
- Ensure that Relief Zone 1 forms the middle layer of the sole!
- Insert the individually adapted sole system into the shoe.
- After putting on a bandage or a sock, the shoe can now be fitted to the patient (Fig. 7).
- The hook and loop fasteners can be shortened as required.

Typical application: patient heavier than 90 kg

Cutting with pressure redistribution zones covering two puzzle parts:
If cutting the relief zones should cover two puzzle parts, it is advisable to anchor the two parts with adhesive tape before cutting at the join.