Your physician has dispensed to you a stabilizing ankle foot orthosis. The design of the rocker sole warrants a physiologically correct gait by means of the »rolling-off« movement. Walking with the FXPro® Walker / AirTraveler™ Walker is more difficult than walking with conventional footwear and its correct application is essential for the Walker's functionality – please read the following precautions.

## Caution

- > Do not attempt to drive while wearing the walker.
- > Use extreme caution when climbing stairs, stepping up onto curbs and walking on uneven surfaces.
- > Walk more slowly than usual. You are not wearing a normal street shoe and more clearance is required for the wedge sole.

The application of this aid should be supervised by a physician or healthcare specialist. In case of pain, swelling or skin irritation, please call a doctor immediately!

### AirTraveler™ Walker only

Avoid excessive compression of the air bladders!

## Indications

- > Medication of edema
- > Acute care in case of ankle joint sprain
- > Injuries of the tissue surface of the lower leg
- > Stress fractures
- > Foot and ankle joint fractures
- > Metatarsal fractures
- > Achilles tendon injuries
- > Injuries of the forefoot, midfoot and hindfoot

### Contraindication (AirTraveler™ Walker only)

The AirTraveler™ Walker is not recommended in case of more than usual foot and ankle joint sensitivity.

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Made in China

Shanghai 201108, China  
No. 128 Lane 3458 Hu Min Road  
DARCO Int. Trading (Shanghai) LTD


Manufactured by:

www.darcointernational.com  
Fax 304.522.0037  
Phone 304.522.4883  
Toll-free 800.999.8866  
Huntington, WV 25701, USA  
810 Memorial Blvd.  
DARCO International, Inc.


info@darco.de, www.darco.de  
Fax +49 (0) 88 07 92 28-22  
Telefon +49 (0) 88 07 92 28-0  
82399 Raisting, Germany  
Gewerbegebiet 18  
DARCO (Europe) GmbH

EC REP

## FXPro® Walker

## AirTraveler™ Walker

1.

**FXPro® Walker**  
(stabilizing ankle foot orthosis)  
**AirTraveler™ Walker**  
(stabilizing ankle foot orthosis)

1. Open the buckles and remove softgood from the Walker.
2. Place the softgood around the leg and close it with the buckles.
3. For placing the foot correctly into the Walker take both hands and spread open the medial and lateral spatula stirrups.
4. After entering the Walker, remove the white security straps at the inside of the stirrups (AirTraveler™ Walker) and subsequently fix the stirrups at the softgood. Then close the buckles (start with the lower buckles). Ensure that the leg fits comfortably into the Walker.

**AirTraveler™ Walker only**

5. Inflate compression bladders of the softgood with air. Use the transparent side of the pump (picture A). The air is distributed evenly. The compression to be realized should not permit free movement of the foot or ankle joint.

The softgood of the Walker is intended to adequately safeguard the foot and ankle joint. In case of tingling, burning or numbness reduce the compression by releasing air from the compression bladders. Utilize the black side of the pump (picture A). If the pain is continuing, interrupt the treatment and consult a doctor.

**Inserting the foam pads**

The four hexagonal pads (picture B) provide immobilization and upholstery of the ankle joint. Place pads where rubbing should be prevented, or where extra space should be filled.

**Cleaning**

For cleaning the shell and the softgood use lukewarm water and possibly mild detergent. Immediately rinse and dry the softgood. Let it air dry. Do not use dryers with warm or hot air.

**Content**



Stabilizing orthosis  
Pump (AirTraveler™ Walker only)  
4 x foam pads

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