

# DARCO Arctic Air® cold compression therapy system

quick and easy to use

highest mobility



detachable cold packs



Arctic Air® Shoulder



Arctic Air® Wrist



Arctic Air® Knee



Arctic Air® Ankle

#### Indications

- > pain therapy
- > rheumatic joint pains
- > chronic arthritis
- > post-surgery recovery
- > joint, ligament and muscle injuries
- > inflammation
- > first aid and rehabilitation



# DARCO






- > immediate pain relief
- > easy to use
- > prompt and simple application
- > preserves movability
- > washable and reusable
- > breathable high-quality tissue
- > incl. hand pump for quick compression and Luer-Lock closure system

## Arctic Air® cold compression therapy system

Arctic Air® has been developed in close cooperation with physiotherapists, physicians and users. The system is designed for pain therapy, for the use after sports injuries and postoperative care. Tests confirm a high cooling efficiency, protective support against influences from outside as well as the preservation of maximum mobility compared to standard products. Hand pumps (with Luer-Lock) are able to generate individual compression. In case of acute care this helps to lower »hemorrhaging« as well as to reduce the length of therapy. The Arctic Air® bandages by DARCO are made exclusively of high-quality materials. They are easy to clean and meet maximum quality standards.

**Attention:** Cold packs to cool in the refrigerator (3 – 8 °C). Through the combination of pressure and cold, the temperature control is recommended in the fridge in any case, it can cause frostbite. In case of sensitivity to cold we recommend to place a cloth between your skin and the bandage. If you have any questions please contact your physician or physiotherapist.



The Arctic Air® bandage protects the injured part of the body

The lower skin layers are able to absorb coldness with immediate effect

Swelling decreases which reduces pain efficiently

### The Arctic Air® technology

The cold pack adjusts perfectly to the body surface. Cold is able to reach lower layers of the skin and consequently reduces the swelling. Pain will recede efficiently right away. When inflating the air cushion compression is generated which increases the effectiveness of therapy. The bandage made from durable fiber offers additional protection.

#### Arctic Air® Cold Packs

The cold packs are detachable and available separately.



#### Arctic Air® Cold Pack | Extra Cold Pack

The extra cold pack guarantees prompt and flexible use of the Arctic Air® cold compression system. The cold pack can be unlimited re-used.

#### Arctic Air® Instant Cold Pack

The instant cold pack will become cold in seconds. It can be activated anywhere, anytime. The instant cold pack can be used once.



**No fridge needed!**

Product	Item-No.	Size
Arctic Air® Ankle	KKT-A	Universal
Arctic Air® Knee	KKT-K	Universal
Arctic Air® Wrist	KKT-W	Universal
Arctic Air® Shoulder	KKT-S	Universal
Arctic Air® Cold Pack Ankle	KKT-AR	Universal
Arctic Air® Cold Pack Knee	KKT-KR	Universal
Arctic Air® Cold Pack Wrist	KKT-WR	Universal
Arctic Air® Cold Pack Shoulder	KKT-SR	Universal
Arctic Air® Instant Cold Pack Ankle	KKT-AI	Universal
Arctic Air® Instant Cold Pack Knee	KKT-KI	Universal
Arctic Air® Instant Cold Pack Wrist	KKT-WI	Universal
Arctic Air® Instant Cold Pack Shoulder	KKT-SI	Universal