

Heelift® AFO Standard and Smooth Ankle Foot Orthosis

Offers the benefits of an ankle foot orthosis device. It provides control of the ankle's position and motion. In addition it allows for limited weight bearing.




The combination of added protection and the ability to customize the product separates Heelift® AFO from other products in the industry.

- > dorsiflexion ankle support for paralyzed legs
- > treatment and prevention of plantar flexion contracture
- > treatment and prevention of foot drop deformity
- > comfortable night splint
- > treatment of plantar fasciitis
- > treatment of achilles tendonitis



DARCO

Heelift® AFO, like all Heelift® Boots, allows for patient-specific customization. The semi-rigid polypropylene brace is designed to provide the best possible control of the ankle's position and motion, while also allowing for limited weight bearing. The Heelift® AFO has a non-slip traction sole that prevents slippage. This feature is ideal for transfers from bed, wheelchairs, commode or physical therapy. Assistance is recommended.



Semi-rigid polypropylene brace

While many orthotics are rigid and put patients at risk, the soft foam of Heelift® AFO, together with the semi-rigid polypropylene brace, provides the best possible control of the ankle's position and motion while also allowing for limited weight bearing. The low friction tricot backing allows better mobility for patients while in bed.

Customizable foam pads

Like Heelift® Suspension Boot, Heelift® AFO allows for customization via a 5 centimeter, beveled fixed pad. The pad elevates the heel and also protects the Achilles tendon.

Additional spare pad

An additional spare pad with an adhesive backing is also included. The spare pad can be used to control hip rotation and malleolar decubitus.

Like other Heelift® products, Heelift® AFO can be washed safely in a net or mesh and can be autoclaved for sterility.



Two types of foam

Heelift® AFO is available with a smooth foam interior (ideal for edematous legs) or a convoluted foam interior (a slightly cooler option).

Item-No.	Colour	Size
Heelift® AFO Standard (convoluted foam):		
DM-HL-AFO	Beige	Universal
Heelift® AFO Smooth (smooth foam):		
DM-HLS-AFO	Beige	Universal

right and left fit