

# An ounce of Prevention and a pound of Cure

As the occurrence of pressure ulcers becomes more prevalent, the pressure is on to find a cost-effective method to prevent them. In the U.S. alone, 1.7 million patients develop pressure ulcers at costs nearing \$3.6 billion annually! Factors in the rise of pressure ulcers presumably will increase as the population ages and obesity-related diseases such as diabetes continue to flourish.

## Heelift® Suspension Boots

- Choice of Original (convoluted foam) or Smooth designs
- Friction-free tricot backing for improved patient mobility
- Extra pad to control hip rotation or foot drop or to provide added elevation
- Polyurethane stiffener to prevent boot from buckling in use
- One size fits all



For non-edematous patients.



For edematous legs or sensitive skin.

Heelift® Original and Smooth Patent No. 5449339

## Elbowlift® Suspension Pad

IMPROVED!

- Cushions and elevates elbow to minimize surface contact
- Smoother, non-friction backing for improved patient movement
- Less bulky for optimum patient comfort
- Stays in place



Patent No. 6216268

**When You Use Heelift®  
and Elbowlift®  
The Pressure's Off!**

For information and pricing, contact our European licensed manufacturer:

**DARCO**

**Darco (Europe) GmbH**  
Gewerbegebiet 18  
82399 Raisting  
Tel: +49-8807-9228-0  
Fax: +49-8807-9228-22  
info@darco-europe.com  
www.darco-europe.com