


The original PegAssist™ System by DARCO for effective off-loading.



- > removable pegs for effective off-loading
- > shock-absorbing effect
- > available for the shoe models shown on reverse side
- > ambilateral use

DARCO



Pegs can be removed from specific areas of the insole in order to redistribute pressure. The cardboard cover then stabilizes the off-loading area.



MedicalSurgical™-ST



OrthoWedge™/ light



SlimLine™



HeelWedge™



A.P.B.™



Walker FXPro™

The 18 mm thick insole is layered with Poron, Plastizote and EVA to provide the perfect combination of shock absorption and moldability.

After removing the pegs from specific areas of the insole, the adhesive-backed cardboard cover is then applied to the complete surface of the insole. The cardboard cover acts to stabilize the remaining pegs.

The PegAssist™ System is available for all products listed here.

MedicalSurgical™-ST

Item-No.	Size (D/EU)	
Women:		
PAS-W1	Small	34,0 – 36,5
PAS-W2	Medium	37,0 – 38,5
PAS-W3	Large	39,0 – 41,0
Men:		
PAS-M1	Small	38,0 – 41,0
PAS-M2	Medium	41,5 – 43,0
PAS-M3	Large	43,5 – 45,0
PAS-M4	X-Large	45,5 – 47,0

All other models

Item-No.	Size (D/EU)	
PAS-0	X-Small	33,5 – 37,0
PAS-1	Small	37,5 – 39,5
PAS-2	Medium	40,0 – 41,5
PAS-3	Large	42,0 – 44,0
PAS-4	X-Large	45,5 – 47,0