

WCS™ Light Closed Wound Care Shoe System



New at DARCO!



- > for treatment of the diabetic foot syndrome
- > for healing of open wounds and minor ulcerations
- > for pressure relief
- > wide toe box for any remaining bandages



DARCO

Pressure points and foot ulcers occur most frequently with diabetes mellitus and in many ways present a great challenge in healing and off-loading sensitive feet.



- > seamless toe area
- > elastic and breathable upper
- > extra space for light bandages in toe area
- > 4 removable insoles in different colours and shore grades; the individual customization of the insoles enables a selective creation of plantar off-loading zones
- > comfortable fitting provided by high heel cup
- > decent look

WCS™ Light Closed Wound Care Shoe System

The diabetic patient requires specifically at an early stage of diabetic foot syndrome with minor ulceration a comfortable and reliable shoe ensuring effective pressure relief for indoor and outdoor use.

The new WCS™ Light of DARCO offers the proven four-layer-insole-system with individual customization for off-loading included in a closed shoe. The ultra soft and breathable upper provides maximum comfort and leaves extra space for bandages.



The WCS™ Light is available pairwise in seven sizes.

Item-No.	Colour	Size	
WCSL0	Black	XS	35.0 – 36.0
WCSL1	Black	S	37.0 – 38.0
WCSL7	Black	S/M	39.0 – 40.0
WCSL2	Black	M	41.0 – 42.0
WCSL3	Black	L	43.0 – 44.0
WCSL4	Black	XL	45.0 – 46.0
WCSL5	Black	XXL	47.0 – 48.0