

# Air Pump Walker®

## Lower Leg Foot Orthosis

Immobilisation in a given position with simultaneous mobilisation

New



Stay mobile for sure!



DARCO

# Features

## Pump system

### COMFORTABLE, EASY AND FAST TO USE

- > Soft lining with integrated air pump
- > Reduction of pain and swelling thanks to the air chamber system
- > Perfect fit shape
  - Safe stabilisation of the ankle



## Construction

### ERGONOMIC DESIGN FOR HIGH WEARING COMFORT

- > Stable, lightweight and breathable
- > Individually adjustable hook and loop straps
- > Flexible calf cuff
- > Ventral cover
- > Hexagonal pads for immobilisation and additional cushioning of the ankle
- > Individual adaptability
- > Adjustable and removable toe guard



## Outsole

### MOBILITY IN DAILY LIFE

- > Slip-resistant and shock absorbing TPR (Thermoplastic Rubber) outer sole for high traction
- > Biomechanically optimised ball rocker
- > Physiologically correct, safe gait pattern



# Three versions

Right = Left 

Low weight 



## Air Pump Walker®

### INDICATIONS

- > Pre- and postoperative / post-traumatic immobilisation in case of damage to the bone /capsule and ligament structures of the foot /ankle
- > Soft tissue injury
- > Stable foot and /or ankle fractures
- > Severe sprains grade II and III



## Air Pump Walker® (Short version)

### INDICATIONS

- > Pre- and postoperative / post-traumatic immobilisation in case of damage to the capsule and ligament structures of the foot /ankle
- > Soft tissue injury
- > Stable foot fractures
- > Sprains grade II



## Air Pump Walker® Achill

### INDICATIONS

- > Achilles tendon ruptures

## Article List



**Air Pump Walker®**  
Lower Leg Foot Orthosis



**Air Pump Walker®**  
Lower Leg Foot Orthosis



**Air Pump Walker® Achill**  
Achilles walker (incl. Achilles Wedges)

SIZE	SHAFT HEIGHT	ART. NO.	SIZE	SHAFT HEIGHT	ART. NO.	SIZE	SHAFT HEIGHT	ART. NO.			
Colour: Grey-Blue			Colour: Grey-Blue <b>Short version</b>			Colour: Grey-Blue					
XS	34.0 – 36.0	34.4 cm	APW0	XS	34.0 – 36.0	21.7 cm	APWS0	XS	34.0 – 36.0	34.4 cm	APWA0
S	37.0 – 39.0	38.0 cm	APW1	S	37.0 – 39.0	23.0 cm	APWS1	S	37.0 – 39.0	38.0 cm	APWA1
M	40.0 – 42.0	41.8 cm	APW2	M	40.0 – 42.0	25.1 cm	APWS2	M	40.0 – 42.0	41.8 cm	APWA2
L	43.0 – 45.0	44.7 cm	APW3	L	43.0 – 45.0	26.4 cm	APWS3	L	43.0 – 45.0	44.7 cm	APWA3
XL	46.0 – 48.0	47.0 cm	APW4	XL	46.0 – 48.0	30.0 cm	APWS4	XL	46.0 – 48.0	47.0 cm	APWA4

## Inclusive

**Liner Sock**  
Stocking

SIZE	ART. NO.
Colour: White Universal (1 pair)	LSOCK

## Combination Products



**APPAW (Air Pump Walker®)**  
Achilles Wedges  
Heel Wedges for Walkers


SIZE	ART. NO.
Colour: Black Universal	APPAW



**TwinShoe**  
Height Compensation Shoe

SIZE	ART. NO.
Farbe: Schwarz	
XS 33.5 – 37.0	TS0B
S 37.5 – 39.5	TS1B
M 40.0 – 41.5	TS2B
L 42.0 – 44.0	TS3B
XL 44.5 – 47.0	TS4B

## Sizes and Measurements



All measurements are given in cm and correspond to the size of the insole. The specified shoe sizes are guidelines. The choice of shoe size depends largely on the volume of the respective dressing.