

Body Armor[®] Pro Term

Foot Stump Orthosis

New
Version



Now
25%
lighter*



NEW
PowerPods[®]
technology

The perfect
solution
**for forefoot
amputations.**

Interim care following amputations.

The Body Armor® Pro Term Foot Stump Orthosis adapts flexibly to stump fluctuations during the interim phase, thanks to a soft lining with integrated pump. The ankle-crossing shaft in a 90° position stabilises the upper and lower ankle, while the ventral support strap and the variable closing system allow for even compression.



NEW Closing system with guiding tabs

Easier closure through better force transmission

Features

Stability and protection

- Shell construction in 90° position
 - Stabilisation of the ankle
- Secure closing system
 - Firm fit and safe stump guidance
- High design with wide tread outsole surface
 - Greater stability and secure fit

Biomechanics and relief

- Soft foot bed sandwich insole
 - Pressure distribution for decubitus prophylaxis
- Orthograde build-up
 - Natural posture and optimised biomechanics

Customised fit and comfort

- 25 % lighter (approx. 220 g) than the previous model
 - Lower load
- Extended stump length (+3 cm)
 - More flexibility in care
- Customisable PowerPods® outsole
 - Adjustment to individual needs
- Flexible ventral support flap
 - Ergonomic fit, even pressure distribution
- Guiding tabs
 - Easier closing of the straps

Practical advantages

- Inspection and ventilation window
 - Better ventilation, improved limb control
- Wide opening
 - Easy to put on / take off
- Universal fit for right and left
- Colour: Grey-Blue

Indications

Postoperatively for partial foot amputations (e.g. Chopart or Lisfranc amputations) due to

- Diabetic foot syndrome
- PAOD (peripheral arterial occlusive disease)
- Trauma



Wide opening for easy foot insertion



Perfect fit through air regulation

Sizes + Measurements

Optimized fit: Please note the revised dimensions.

Size	Length Stump	Insole in cm			Art.-Nr.
		B	C	E	
M	19.5 cm	13.0	9.5	20.0	BAPT2
L	21.0 cm	14.0	10.5	22.0	BAPT3
New XL	22.5 cm	14.5	11.5	23.6	BAPT4



Measurements and sizes are approximate and depend on the volume of the dressing applied.

Compatible product

WCS®

Wound Care Shoe



Right = Left



For diabetics



Low weight



Time saving



DARCO (Europe) GmbH

Gewerbegebiet 18 | 82399 Raisting | Germany | Tel. +49 8807 9228-0 | info@darco.de | www.darco.de

