MedSurg Postoperative Shoe for Stabilizing and Off-loading





The DARCO MedSurg Shoe is one of the most versatile postoperative shoes on the market today. The original insole can be customized to better meet the needs of individual patients.






The MedSurg Shoe with its famous semi-rigid sole is built around our patented MetaShank™ forefoot support, that provides more rigid control and better protection after osteotomies.

Our breathable, mesh upper conforms to bony abnormalities of the foot and has a reinforced heel counter and pad ded heel collar to stabilize the rear foot and reduce lippage. Extra-long straps accommodate even the bulkiest dressings. Due to the classification offered in men's and women's sizes the shoe fits at best.

Indications

> used postoperatively for stability and off-loading, e. g. after foot osteotomy procedures





Features and Benefits

- > Square Toe Design acts as a bumper to provide extra protection and facilitates better left/right fit
- > MetaShank™ protection provides more rigid control under the metatarsal heads
- > Stapless Closure eliminates buckle pressure on the forefoot and can expand to accommodate dressing up to 17 inches in circumference
- > Rocker Sole reduces plantar pressure on the heel and forefoot by over 25%
- > Fully compatible with DARCO's PegAssist® Off-loading Insole or Puzzle Insole® Modular Insole



SIZE

Women | Colour: Black

WS	34.0 - 36.5	MSW1B-ST
WM	37.0 - 38.5	MSW2B-ST
WI	39 0 - 41 0	MSW3R-ST

Men | Colour: Black

MS	39.0 - 41.0	MSM1B-ST
MM	41.5 - 43.0	MSM2B-ST
ML	43.5 – 45.0	MSM3B-ST
MXL	45.5 - 47.0	MSM4B-ST

Children | Colour: Black up to 32.0



ART. NO.

Optional







Puzzle Insole®



PegContour® Insole



PegAssist® Insole



Relief Contour Insole

