


OrthoWedge Pro with EVA Outsole

Forefoot Off-Loading Shoe

Forefoot off-loading after surgery to correct hallux valgus or hammer toes, and for the healing of forefoot ulcerations

20%
lighter
EVA outsole *



For faster and better regeneration

- > Light weight
- > EVA outsole: Abrasion-resistant with excellent damping properties
- > Easy replacement of the insole

* Comparison with industry-standard post-op shoes with TPR outsole



DARCO

The OrthoWedge Pro is the optimised version of the OrthoWedge



Features

- > EVA outsole → Abrasion-resistant with excellent damping properties
- > Widened footprint
- > Wide hook and loop strap with medial and lateral strap loops
- > Soft, seamless liner → High wearing comfort
- > Pressure redistribution to the midfoot or rearfoot
- > Specific forefoot off-loading
- > Square Toe: Generous angular shape → Toe protection
- > Wide opening and adjustable straps → Comfortable entry
- > Easy modification to the right or left foot

Indications

- > Postoperative use after surgery to correct hallux valgus, hammer / claw toes
- > Healing of forefoot ulcerations

Offloading in the treatment of a diabetic foot ulcer is only indicated when other necessary offloading measures (e.g. a total contact cast or a comparable orthosis) are definitively not possible.

Sizes and Measurements

	Size	Shoe size	Insole in cm			Art. No.
			A	B	C	
	XS	33.5 – 37.0	21.5	8.5	6.5	OW0B-STL
	S	37.5 – 39.5	23.5	9.5	7.0	OW1B-STL
	M	40.0 – 41.5	25.5	10.0	7.5	OW2B-STL
	L	42.0 – 44.0	27.5	10.5	8.5	OW3B-STL
	XL	44.5 – 47.0	29.5	11.5	9.0	OW4B-STL

Compatible with

