Relief Dual® Plus

Off-loading shoe with high ankle support

PLUS even more applications for postoperative management and wound treatment

ap to size MXXL (47.5 – 50.0)

64 %
higher than the Relief Dual®

> The Relief Dual® Plus easily accommodates bulky dressings

- > For postoperative treatment and for wound therapy
- > Relief and restriction of ankle movement thanks to shoe base reinforcement
- > Innovative dual sole technology: Abrasion-resistant, light, shock-absorbing and anti-slip
- > Can be combined individually with DARCO modular insole systems






Relief Dual® Plus

Off-loading shoe with a high shaft and adjustable heel cap













Proven quality, innovative concepts!Back-shifted pivot point for less pressure on the forefoot





Neutral

Relief Dual®/Plus

Flat off-loading shoe with a high shaft and a closure strap that goes around the ankle

There is virtually no movement of the ankle during each step due to continuous, rigid reinforcement of the base of the shoe which generates greater off-loading of the forefoot. Furthermore, this achieves improved pressure redistribution.

INDICATIONS

- > Postoperative
 - For stabilisation and offloading, e.g. foot osteotomies
 - · For removal of Haglund exostosis
 - · After surgery for fasciitis plantaris / heel spur
- > For wound therapy

FEATURES

- > Shaft: 64% higher than in the Relief Dual®
- > High shaft → Orthograde rearfoot stability
- > Individually adjustable heel cap → Optimal support + maximum stability
- > Closure strap goes around the ankle → Additional support
- > Maximum opening width → Space for bulky dressings
- > Innovative dual-sole technology → Abrasion-resistant, shock-absorbing, anti-slip, with continuous rigid shoe base reinforcement
- > Reduced dorsal extension momenta → High forefoot off-loading
- > Less ankle movement + improved pressure redistribution after surgery
- > Negative heel through use of the insoles Puzzle Insole®, Relief Contour Insole
- > Suitable for diabetes-adapted footbed
- > Wide hook and loop strap with medial + lateral strap loops
- > Back-shifted ball rocker → Less loading in the forefoot
- > Square Toe: Generous angular shape → Toe protection

Benefits at a glance



- > High shaft for orthograde rearfoot stability
- > Flexible and adjustable heel cap
- > More space, no pressure points
- > Soft lining



- > Closure strap goes around the ankle
- > Strap loops on both sides
- > Hook and loop fastener can be cut to desired length



> Latest rear entry design due to maximal opening of the shoe



> Seamless and larger comfort area after, e.g., hallux valgus surgery



> Built up in layers for improved pressure redistribution through the long-soled, rigid inner sole



- > Space for inserts through the removable sole (also diabetes-compliant foot bed)
- > Low profile height no height adjustment necessary



- > Slip-resistant rubber sole
- > Improved shock absorption particularly in the heel and ball area



- > Toe cap optionally available
- > Impact protection for the toes
- > Protects toes from inclement weather and dirt



Relief Dual® Plus, PegContour® Insole Uni and Toe Cap

> Available p to size MXXL (47,5 – 50,0)



Relief Contour Insole

Arch Cushion Insole



For diabetics



Time saving







Pressure redistribution of up to 90% due to nearly full contact

INDICATIONS

- > After hallux valgus and small toe surgery (Austin, Akin, Scarf, Weil, Hohmann etc.)
- > Talipes valgus, flat foot, splay foot
- > Angio-neuropathic foot changes in combination with mild foot deformities

FEATURES

- > Arch cushion insole (long-soled)
- > Anatomically adapted contours
 - → Optimal pressure redistribution
- > Contact between foot and sole over almost the entire surface
 - → Pressure redistribution up to 90%
- > Dorsal extension: Negative heel of 2.4° \rightarrow Forefoot off-loading
- > Medial longitudinal arch and proximal transverse arch
 - → Ideal load redistribution
- > Customisable PU material with high shape recovery
- > Soft lengthwise and crosswise sole covering
 - → High wearing comfort
- > Square Toe: Generous angular shape
 - → Toe protection





PegAssist® Insole

Peg Off-Loading Insole







For diabetics










Relief Dual®/Plus



PegAssist® Insole (MTP II + III omitted) in the Relief Dual®/Plus









INDICATIONS

For selective off-loading for:

- > Plantar lesions
- > Plantar ulcerations
- > Pressure points

FEATURES

- > Insert system with individually removable pegs
 - → Customisable fitting
- > Plantar pressure redistribution
- > Targeted, punctual pressure relief
- > Pressure reduced by up to 50%
- > Shock-inhibiting
- > Material: 18 mm PE foam layer + Poron® + EVA
- > Square Toe: Generous angular shape
 - → Toe protection
- > Marker pen (included)



PegContour® Insole Uni

Universal Peg Arch Cushion Insole



For diabetics











Due to the symmetry, the insole is wearable either on the right or left.





(MTP II + III omitted) in the Relief Dual®/Plus







Simultaneous targeted offloading and pressure redistribution over the entire surface

INDICATIONS

- > After hallux valgus and small toe surgery (Austin, Akin, Scarf, Weil, Hohmann etc.)
- > Postoperative use and treatment of the heel
- > Angioneuropathic foot changes in combination with mild foot deformities
- > Healing of ulcerations in the heel area (rearfoot ulcers)
- > Fasciitis plantaris / heel spur
- > Heel wounds

For selective off-loading for:

- > Plantar lesions
- > Plantar ulcerations / ulcer prevention
- > Pressure points

FEATURES

- > Shock-inhibiting arch cushion insole (long-soled) with individually removable pegs → Customisable fitting
- > Plantar pressure redistribution + targeted, selective pressure relief
- > Pressure reduced by up to 45%
- > Customisable PU material with high shape recovery
- > Soft sole covering made of lengthwise and crosswise elastic polyester
 - → High wearing comfort
- > Symmetrical sole
 - → Universal fit for right + left foot
- > For lesions in the arch of the foot
- > Anatomically adapted contours
 - → Optimal pressure redistribution



The new version have several advantages:

- + Modified production method
- → No EVA insole floor / backing plate required
- + Easier and quicker customisation
- + Optimised PU material composition
 - → Better resilience
- + New insole lining and colour design
- + Made in Europe

Puzzle Insole®

Modular Insole

















Indication-specific pressure relief through a modular insole with interchangeable elements

INDICATIONS

Relief Dual®/Plus

1 Small Toe Module (narrow)

> Tailor's bunion (digitus superductus)

(big toe module changed) in the Relief Dual®/Plus

> Jones fractures (puzzle module extends to Lisfranc joint row)

2 Big Toe Module (wide)

- > Metatarsophalangeal joint arthrodesis
- > Austin, Scarf and Akin osteotomies
- > Closing and opening wedge
- > Transverse and pendulum osteotomies
- > Lapidus arthrodesis (TMT1 arthrodesis)
- Youngswick osteotomies (puzzle module extends to cuneiform I)

3 Middle Module

- > Weil and/or Helal osteotomies I IV
- > Morton's neuroma
- > Hammer toe corrections
- > PIP/DIP arthrodesis
- > Metatarsal bone fractures II IV
- > Shortening osteotomies of the metatarsal bones II - IV (puzzle module extends to proximal end of the metatarsal bones II - IV)







Turn two into one to customize!

FEATURES

- > Set of 2 soles (blue + black) with 4 elements each → Individual combination options
- > Different Shore hardnesses: Blue 09 | Black 30
- > Different composition of the sole modules
 - → Indication-dependent pressure redistribution
- > Dorsal extension: Negative heel of 2.4°
 - → Forefoot off-loading
- > Plantar pressure redistribution

Shore hardness

- > Shock-inhibiting
- > Material: EVA foam
- > Square Toe: Generous angular shape
 - → Toe protection

Article List





PegContour® Insole Uni



SIZE		ART. NO.	SIZE		ART. NO.	SIZE		ART. NO.	
Woman Colour grey			Woman			Woman			
WS	34.0 – 36.5	RDP-W1	WS	34.0 – 36.5	PCIU-W1	WS	34.0 – 36.5	TCDW1	
WM	37.0 – 38.5	RDP-W2	WM	37.0 - 38.5	PCIU-W2	WM	37.0 – 38.5	TCDW2	
WL	39.0 – 41.0	RDP-W3	WL	39.0 – 41.0	PCIU-W3	WL	39.0 – 41.0	TCDW3	
Men Colour grey			Men			Men			
MS	39.0 – 41.0	RDP-M1	MS	39.0 – 41.0	PCIU-M1	MS	39.0 – 41.0	TCDM1	
MM	41.5 – 43.0	RDP-M2	MM	41.5 – 43.0	PCIU-M2	MM	41.5 – 43.0	TCDM2	
ML	43.5 – 45.0	RDP-M3	ML	43.5 – 45.0	PCIU-M3	ML	43.5 – 45.0	TCDM3	
MXL	45.5 – 47.0	RDP-M4	MXL	45.5 – 47.0	PCIU-M4	MXL	45.5 – 47.0	TCDM4	
741742									
MXXL	47.5 – 50.0	RDP-M5	MXXL	47.5 – 50.0	PCIU-M5	MXXL	47.5 – 50.0	TCDM5	



Relief Contour Insole



PegAssist® Insole



Puzzle Insole®

SIZE		ART. NO.	ART. NO.	SIZE		ART. NO.	SIZE		ART. NO.
Women		Right	Left	Wome	Women		Wome	Women	
WS	34.0 – 36.5	RC-W1-R	RC-W1-L	WS	34.0 – 36.5	PAS-W1	WS	34.0 – 36.5	PI-WS
WM	37.0 – 38.5	RC-W2-R	RC-W2-L	WM	37.0 – 38.5	PAS-W2	WM	37.0 – 38.5	PI-WM
WL	39.0 – 41.0	RC-W3-R	RC-W3-L	WL	39.0 – 41.0	PAS-W3	WL	39.0 – 41.0	PI-WL
Men		Right	Left	Men			Men		
MS	39.0 – 41.0	RC-M1-R	RC-M1-L	MS	39.0 - 41.0	PAS-M1	MS	39.0 – 41.0	PI-MS
MM	41.5 – 43.0	RC-M2-R	RC-M2-L	MM	41.5 - 43.0	PAS-M2	MM	41.5 - 43.0	PI-MM
ML	43.5 – 45.0	RC-M3-R	RC-M3-L	ML	43.5 – 45.0	PAS-M3	ML	43.5 – 45.0	PI-ML
MXL	45.5 – 47.0	RC-M4-R	RC-M4-L	MXL	45.5 – 47.0	PAS-M4	MXL	45.5 – 47.0	PI-MXL

Sizes and Measurements

















