

TAS[®] (Toe Alignment Splint)

Toe Alignment Splint



Elastic heel strap

Elastic metatarsal strap



High wearing comfort, even under socks

Universal size



Right = Left



Low weight



Latex free



DARCO

The TAS® Toe Alignment Splint is a bandage for alignment. It is used postoperatively instead of taping and is an effective and affordable way to fix the toes in the correct position after surgery / correction and to maintain the desired toe position postoperatively during the healing process. The elastic heel strap at the back of the foot reduces slippage of the midfoot bandage and the individual loops keep the toes in the desired position. The soft material provides a high wearing comfort. The TAS® is also an optimal care after minimally invasive surgeries.

Indications

Postoperative following:

- 1) Hallux valgus surgery
- 2) Hammer toe surgery
- 3) Tailor's bunion surgery



Eigenschaften

- > Postoperative bandage for alignment used instead of taping
- > Freely positionable fixation + alignment of operated / corrected toes
- > Also for use after minimally invasive procedures
- > Soft, skin-friendly material → Pleasant wearing comfort
- > Stabilising element thermoplastic deformable
- > Elastic, 34 cm long metatarsal strap → Care of wide feet possible
- > Elastic heel strap on the rearfoot
- > Can be worn with post-op and healing shoes
- > Universal fit for right and left



Adduction



Abduction



Dorsal extension



Plantar flexion



SIZE ART. NO.



Colour: White
Universal TAS-22