



DARCO

Ubiq[®]

NEW

The everyday hero.
Style included.

Therapeutic Comfort Shoe



More lightness in every step.



Loop
for easy donning

Soft
liner

Relief
for the forefoot

Space for
your own insoles

Are you looking for...

... comfortable shoes for
everyday wear?

... styles that work with your
own insoles?

... a look that fits perfectly
at work and in your free time?

Then let us introduce you to the new Ubiq®!

A shoe for anyone who wants to be comfortable,
healthy and stylish on the go. The therapeutic
comfort shoe Ubiq® supports your natural rolling
motion and provides an all-around pleasant walking
experience – all day long.

Whether at work, on a walk, or at home:
Ubiq® moves with you every step of the way.
Functional, modern, and versatile – a true
everyday hero.

Your way. Your mobility. Your Ubiq®.



Indications.

When the Ubiq® provides you with optimum support.

- Mild to moderate foot / toe deformities such as hallux valgus, hammer toes or claw toes
- Hallux limitus
- Diabetes mellitus without PNP / PAOD (risk group 0 according to the German Diabetes Society)
- Diabetes mellitus without PNP / PAOD with foot deformity (risk group 1 according to the German Diabetes Association)

Features.

What makes the Ubiq® special.



Outsole

- Back-shifted ball rocker
 - Forefoot off-loading
- Integrated, large, abrasion-resistant rubber coating (4 mm)
 - Non-slip
- Shock-absorbing, force-absorbing sole construction
- Continuous and reinforced rigid shoe base
 - Forefoot off-loading
- Reduced footprint
 - Dynamic, fast rolling motion
- Lower heel-to-toe drop (7 mm)
 - Natural foot position



Shoe Last Ubiq®

- Foot-friendly last shape
- Ideal ball width
- Tailored shape
- Toe box height 23.3 – 27.6 mm (depending on size)



Insole

- Anatomically moulded footbed (6 mm)
 - Moderate heel guide
 - Slight lengthwise and crosswise arch support
 - PU sole with smooth microfibre cover for good cushioning + easy foot insertion
- Replaceable
- Sufficient space (6 mm) for customised insoles




Lacing

- Elasticated laces with quick-release fastener
 - Easy to put on and take off
- Non-elastic, flat laces (spare set enclosed)
 - Better support



Shaft

- Air-permeable
- Reinforced toe
- Anatomical ankle cutout
 - Ideal rearfoot fit
- Optimally adapted heel cup
 - Substantial and cushioned heel support
- Loop on the heel
 - Easy to put on
- Stable upper material
- Soft liner
 - Good padding
 - Protection against pressure and external influences
 - Comfortable to wear



The back-shifted ball rocker relieves pressure on the forefoot

Four styles.

Which Ubiq[®]
are you?



dynamic
shock-absorbing
non-slip

Ubiq[®] Way

Black



Grey



Navy



Lemon



**Available in
shoe sizes 36 – 47**

Find the right
size for you



DARCO (Europe) GmbH

Gewerbegebiet 18 | 82399 Raisting | Germany | Tel. +49 8807 9228-0 | info@darco.de | www.darco.de

